

# Segmentation of Axillary and Supraclavicular Tumoral Lymph Nodes in PET / CT: A Hybrid CNN / Component-Tree Approach


**D. L. Farfan Cabrera, N. Gogin, D. Morland, B. Naegel, D. Papathanassiou, N. Passat**

Université de Reims Champagne Ardenne, CReSTIC, EA 3804, 51097 Reims, France

General Electric Healthcare, Buc, France


Département de Médecine Nucléaire, Institut Godinot, Reims, France

Université de Strasbourg, CNRS, ICube, Strasbourg, France




# Automatic axillary lymph node tumor segmentation in PET/CT

**What are lymph nodes?**



**Patient with no breast cancer**




**Patient with breast cancer**


**lymph nodes are the first  
organs reached by cancer**

# Automatic axillary lymph node tumor segmentation in PET/CT

## What are our contributions?



Patient with no breast cancer



Patient with breast cancer

### 1<sup>st</sup> Contribution:

- Provide a prognostic factor for the staging of breast cancer

**lymph nodes are the first organs reached by cancer**

**Cancer staging is determined by:**


**N:** number of lymph nodes with cancer  
**T:** tumor volume  
**M:** metastasis state

# Automatic axillary lymph node tumor segmentation in PET/CT

## What are our contributions?



PET image



PET image


Brown fat shows up (due to exposure to cold)

### 2<sup>nd</sup> Contribution:

- Help doctors to identify quicker false positive coming from brown fat

**brown fat**

- Brown fat has the same metabolism as tumor cells
- Brown fat generates **false positives** that can be mistaken by tumors




# Automatic lymph node tumor segmentation in PET/CT

## 2. Component-tree extraction

What is a component-tree?

3	3	3	3	3	3	3
3	2	2	2	2	2	3
3	2	1	2	4	2	3
3	2	2	2	2	2	3
3	3	3	3	3	3	3
3	0	0	0	3	3	3
3	3	3	3	3	3	3

Full image



- Each node is a connected component
- Root node contains all the image
- Follows a hierarchical structure according to gray levels
- Each tree level corresponds to a gray level
- Each node belongs to its ascendants nodes

# Automatic lymph node tumor segmentation in PET/CT


## 2. Component-tree extraction



PET image  
Coronal MIP view



Full PET image




# Automatic lymph node tumor segmentation in PET/CT

## 2. Component-tree extraction

For each node  $N_i$ , it is computed:

- $f_1 \rightarrow G(N_i)$ : mean gradient of node contour in PET
- $f_2 \rightarrow H(N_i)$ : mean HU value in CT
- $f_3 \rightarrow S(N_i)$ : standard deviation of H
- $f_4 \rightarrow R(N_i)$ : relative integral volume
- $f_5 \rightarrow L(N_i)$ : position with respect to the lungs



# Automatic lymph node tumor segmentation in PET/CT

## 2. Component-tree extraction

### Goal of these descriptors?

We assign to each **node** 5 descriptors

For each node  $N_i$ , it is computed:

- $f1 \rightarrow G(N_i)$ : mean gradient of node contour  $\rightarrow$  lymph node have **high** contour gradient
- $f2 \rightarrow H(N_i)$ : mean HU value in CT  $\rightarrow$  lymph nodes have **positive** HU value  
brown fat has **negative** HU value
- $f3 \rightarrow S(N_i)$ : standard deviation of H  $\rightarrow$  lymph nodes have **low** standard deviation
- $f4 \rightarrow R(N_i)$ : relative integral volume  $\rightarrow$  lymph nodes have **high** contrast with their neighbourhood
- $f5 \rightarrow L(N_i)$ : position with respect to the lungs  $\rightarrow$  lymph nodes are **outside** the lungs convexhull

# Automatic lymph node tumor segmentation in PET/CT

## 3. Feature map generator

### From feature vectors to feature volumes?


We assign to each **voxel** 5 descriptors

For each voxel  $x$  in the PET image, it is computed:

- $G(x) = \max_i G(N_i)$ 
→ lymph node have
high
contour gradient
- $H(x) = \text{mean}_i H(N_i)$ 
→ lymph nodes have
positive
HU value  
brown fat has
negative
HU value
- $S(x) = \min_i S(N_i)$ 
→ lymph nodes have
low
standard deviation
- $R(x) = \max_i R(N_i)$ 
→ lymph nodes have
high
contrast with their neighbourhood
- $L(x) = \text{or}_i L(N_i)$ 
→ lymph nodes are
outside
the lungs convexhull








# Automatic lymph node tumor segmentation in PET/CT

## 4. CNN and component-tree fusion



# 5. CNN

## 3D U-Net PET-CT for LN tumor segmentation (2 encoders + 1 decoder)



# Automatic lymph node tumor segmentation in PET/CT

## 6. Results

Training: 201 tumors (42 PET/CT exams)  
Validation: 56 tumors  
Patches for test: 63 tumors (10 PET/CT exams)  
Loss Function: 1 - Dice  
Iterations: 1000

Type of CNN: U-NET  
Number of layers: 3  
Resolution:  $1.2\text{mm}^3$   
3D Patch size:  $80\text{ mm}^3$


U-PET → inputs {PET}

U-PET-CT → inputs {PET, CT}

U-PET-CT-all → inputs {PET, CT, all feature maps}

baseline →


proposed  
method →



# Automatic lymph node tumor segmentation in PET/CT

## 6. Results

U-PET → inputs {PET}  
 U-PET-CT → inputs {PET, CT}  
 U-PET-CT-all → inputs {PET, CT, all feature maps}



baseline →  
 proposed method →


3D CNN Model	Voxel-based metrics (Segmentation)			Region-based metrics (Detection)		
	DSC voxel	PPV voxel	SE voxel	DSC region	PPV region	SE region
3D U-Net (PET)	<b>0.832703</b> +- 0.1327	<b>0.827762</b> +- 0.1483	<b>0.82478</b> +- 0.0868	<b>0.798685</b> +- 0.1547	<b>0.757925</b> +- 0.2206	<b>0.890805</b> +- 0.0821
3D U-Net (PET-CT)	<b>0.865896</b> +- 0.08	<b>0.844429</b> +- 0.1034	<b>0.886527</b> +- 0.0465	<b>0.871825</b> +- 0.1312	<b>0.849864</b> +- 0.1893	<b>0.926111</b> +- 0.0628
3D U-Net (PET-CT-all)	<b>0.867845</b> +- 0.1161	<b>0.872781</b> +- 0.1425	<b>0.896473</b> +- 0.0553	<b>0.894624</b> +- 0.1431	<b>0.851361</b> +- 0.1948	<b>0.933135</b> +- 0.0756

# Automatic lymph node tumor segmentation in PET/CT

## 6. Results

Proposed method removes false positives coming from brown fat

Results shown  
on the PET/CT  
fusion images



# Automatic lymph node tumor segmentation in PET/CT

## 6. Results

Proposed method separates visually connected tumors


Results shown  
on the PET/CT  
fusion images

*false positive*


*false positive*

*joined contours*


*separated contours*



U-PET



U-PET-CT  
(baseline)



U-PET-all  
(proposed method)



Thank you for your attention

mail: [diana.lucia.farfan@gmail.com](mailto:diana.lucia.farfan@gmail.com)

page: [linkedin.com/dianafarfan](https://www.linkedin.com/dianafarfan)